

AFM Probe Specifications:

Coating

Reflective Gold

Additional Info

AFM probes of the 160AC series are designed primarily for standard tapping mode AFM imaging in air or vacuum.

The tetrahedral AFM tip is located precisely at the free end of the AFM cantilever. This allows the AFM tip to be positioned accurately over the area of interest on the sample surface.

The uncoated AFM tip offers a sharp AFM tip apex and chemical inertness. The back side gold coating ensures high and stable laser reflectivity in air, vacuum, liquid and aggressive chemical environments.

AFM Tip:

Shape	Height	Setback	Radius	Half Cone Angle
Optimized Positioning	14 μm (12 - 16 μm)*	0 μm	< 7 nm	0° front, 35° back, <9° side

* typical values

AFM Cantilever:

Cantilever	Shape	Force Const.	Res. Freq.	Length	Width	Thickness
Standard tapping mode AFM cantilever	Beam	26 N/m (8 - 57 N/m)*	300 kHz (200 - 400 kHz)*	160 μm (1 - 170 μm)*	40 μm (38 - 42 μm)*	4 μm (3.5 - 4.5 μm)*

* typical values