

AFM Probe Specifications:

Coating

Gold Overall

Additional Info

AFM probes of the 240AC series are designed for tapping mode AFM imaging of soft samples.

The tetrahedral AFM tip is located precisely at the free end of the AFM cantilever. This allows the AFM tip to be positioned accurately over the area of interest on the sample surface.

The gold coated AFM tip is suitable for biological applications, AFM tip functionalization and custom applications. The overall gold coating ensures inertness and electrical conductivity, as well as high and stable laser reflectivity in air, vacuum, liquid and aggressive chemical environments.

AFM Tip:

Shape	Height	Setback	Radius	Half Cone Angle
Optimized Positioning	14 μm (12 - 16 μm)*	0 μm	< 30 nm	0° front, 35° back, <9° side

* typical values

AFM Cantilever:

Cantilever	Shape	Force Const.	Res. Freq.	Length	Width	Thickness
Soft tapping mode AFM cantilever	Beam	2 N/m (0.6 - 3.9 N/m)*	70 kHz (45 - 90 kHz)*	240 μm (1 - 250 μm)*	40 μm (38 - 42 μm)*	2.6 μm (2.1 - 3.1 μm)*

* typical values