

All-In-One-DLC

AFM Tip

SHAPE	HEIGHT	SETBACK	RADIUS	HALF CONE ANGLE
Rotated	17 μm (15 – 19 μm)*	15 μm (10 – 20 μm)*	15 nm	20°–25° along cantilever axis, 25°–30° from side, 10° at the apex

4 AFM Cantilevers

Cantilever A – Contact Mode		Cantilever B – Force Modulation	
Shape	Beam	Shape	Beam
Force Constant	0.2 N/m (0.04 – 0.7 N/m)*	Force Constant	2.7 N/m (0.4 – 10 N/m)*
Resonance Frequency	15 kHz (10 – 20 kHz)*	Resonance Frequency	80 kHz (50 – 110 kHz)*
Length	500 μm (490 – 510 μm)*	Length	200 μm (200 – 220 μm)*
Width	30 μm (25 – 35 μm)*	Width	30 μm (25 – 35 μm)*
Thickness	2.7 μm (1.7 – 3.7 μm)*	Thickness	2.7 μm (1.7 – 3.7 μm)*
Cantilever C – Soft Tapping		Cantilever D – Tapping Mode	
Shape	Beam	Shape	Beam
Force Constant	7.4 N/m (1 – 29 N/m)*	Force Constant	40 N/m (7 – 160 N/m)*
Resonance Frequency	150 kHz (70 – 230 kHz)*	Resonance Frequency	350 kHz (200 – 500 kHz)*
Length	150 μm (140 – 160 μm)*	Length	100 μm (90 – 110 μm)*
Width	30 μm (25 – 35 μm)*	Width	50 μm (45 – 55 μm)*
Thickness	2.7 μm (1.7 – 3.7 μm)*	Thickness	2.7 μm (1.7 – 3.7 μm)*

* typical range

Coating

Diamond–Like–Carbon coating on tip side of the cantilever, 15nm thick;Aluminum coating on detector side of the cantilever, 30 nm thick

Alignment Grooves

none

Additional Info

Versatile silicon AFM probe with 4 different AFM cantilevers on a single AFM holder chip. Rotated, monolithic silicon AFM probe with symmetric tip shape for various applications (contact mode, force modulation mode, tapping mode and soft tapping mode).

High durability and hydrophobicity due to Diamond–Like–Carbon coating on tip side of the AFM cantilever.

The rotated AFM tips allow for more symmetric representation of high sample features. The consistent AFM tip radii ensure good resolution and reproducibility.

The AFM holder chip fits most commercial AFM systems as it is industry standard size.

The main advantage of this product compared to regular, single–cantilever AFM probes is the freedom to choose in the very last moment the right AFM cantilever for each application. You don't need to stock various AFM probe types any more. Nevertheless, this product is not meant as a substitution to comparable single–cantilever AFM probes, because the geometry of each one of the All–In–One cantilevers differs from the geometry of comparable specialized single–cantilever AFM probes.